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ELCON Medical Instruments GmbH
Dr. Karl-Storz-Straße 26
78532 Tuttlingen
Germany

Tel.: +49 7461 9281 - 0
www.elcon-medical.com
info@elcon-medical.de

ELCON
MEDICAL INSTRUMENTS GmbH

English

Processing instructions D

Important notes

Read these reprocessing instructions carefully before each use and keep them easily accessible for the user or the relevant technical personnel.

Carefully read the warnings indicated by this symbol. Improper use of the products may result in serious injury to the patient, users or third parties.

This reprocessing instruction is valid for all article number in Annex I.

Section	Title		Description
10	Intended use and area of application		The devices may only be used for their intended use in medical specialties, by appropriately trained and qualified personnel. The devices are not intended for use on the central circulatory or nervous system. The attending physician or the user is responsible for the selection of the devices for specific applications or operative use, appropriate training and information, and sufficient experience for handling the devices.
20	Precautionary measures and warnings	Treatment of brand new devices	Factory-new devices must have undergone the complete reprocessing process once before being used. Protective caps and protective nets for sharp-edged devices must be completely removed beforehand.
		Load	Overloading of the devices must be avoided. Overloading can lead to bending or breaking and thus to the loss of functionality of the devices.
		Functional impairment	The devices corrode and their function is impaired if they come into contact with aggressive substances. For this reason, it is essential to follow the reprocessing and sterilization instructions.
		Operating conditions	To ensure the safe operation of the aforementioned devices, correct maintenance and care of the devices is essential. In addition, a functional or visual check should be carried out before each use. For this reason, we refer to the relevant sections in these reprocessing instructions.
		Storage	There are no specific requirements for storing the devices before sterilization. Nevertheless, we recommend storing the medical devices in a clean and dry environment .
30	Liability and warranty		The devices may only be used for their intended purpose in the medical fields by appropriately trained and qualified personnel. The attending physician or the user is responsible for the selection of the devices for specific applications or operative use, for appropriate training and information, and for sufficient experience in handling the devices. The manufacturer does not take responsibility for damage resulting from repair or maintenance by unauthorized parties.
40	Obligation to report		All serious incidents related to the devices must be reported to the manufacturer and the competent authority of the member state where the user and/or patient is established.
50		Sterility: Delivery condition	The medical devices are delivered in a non-sterile condition and must be cleaned and disinfected by the user in accordance with the following instructions before the first and each subsequent use.

Section	Title		Description
60	Limitation of reprocessing, disposal		Frequent reprocessing has little effect on the devices. Product life is limited by wear due to intended use and damage to the device. The device may no longer be used, among others, under the following aspects: Corrosion, damage, fractures, cracks, deformation, porosity, functional limitations, devices with unrecognizable or missing labeling. Therefore, the corresponding instructions for functional testing by the reprocessor must be implemented. Please dispose of the devices properly or to a recycling system at the end of the product life. The national regulations and disposal guidelines must be observed!
70	Manual cleaning and wipe disinfection		See the following points
80		Warning notices	Do not immerse the devices in liquids or autoclave them! Tap water to be used must correspond to drinking water quality. If not a one-piece product, the device must be disassembled as far as possible.
90	Manual cleaning		Visible contamination must be wiped off immediately after application. Avoid long waiting times, e.g. overnight or over the weekend (< 6 hours drying time after soiling). If the product is not a one-piece product, it must be disassembled as much as possible. - Remove coarse soiling with a disposable cloth or similar soaked in detergent/disinfectant and dispose of it. - Afterwards, no more contamination should be visible.
100	Wipe disinfection		The special instructions of the disinfectant manufacturer must be observed. Only disinfectants that are suitable for disinfecting the devices being used may be used. (list of the association of applied hygiene, VAH) Wet the device or a lint-free disposable cloth with the disinfectant. Wipe the entire surface of the device. Observe the disinfectant manufacturer's exposure time. Then dispose of the used cloths as residual waste.
110	Drying		Allow the surfaces to air dry.
120	Inspections		The devices must be macroscopically clean after each cleaning, i.e. free of visible dirt. - A stained device (corrosion, discoloration) must be sorted out immediately and subjected to special treatment. - A contaminated device must be reprocessed. - In the event of errors or damage, the device must be sorted out immediately. The following components must be checked with particular attention: - Notches - Blades - Points - If the UDI carrier is illegible, the products must be disposed of.
130		Additional Information	The reprocessor is responsible for ensuring that the reprocessing carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired results. This usually requires validation and routine monitoring of the process and equipment used.
140	Service, repair and return transport	 Service and repair	Do not carry out any repairs or modifications to the devices yourself. This is the sole responsibility of the manufacturer's authorised personnel. Please contact us if you have any complaints, claims or information regarding our devices.
		 Return transport	Defective or non-conforming devices must have gone through the entire reprocessing process before being returned for repair/service.
150	Storage and transport		- Protect against mechanical damage - Store dry and dust-free - Store and transport in safe containers/packaging - Handle with great care, do not throw or drop

ELCON Medical Instruments GmbH
Dr. Karl-Storz-Straße 26
78532 Tuttlingen
Germany

Tel.: +49 7461 9281 - 0
www.elcon-medical.com
info@elcon-medical.de

ELCON
MEDICAL INSTRUMENTS GmbH

Section	Title	Description
160	Material resistance	When selecting cleaning agents and disinfectants, make sure that they do not contain the following ingredients: • Organic, mineral and oxidising acids (minimum permissible pH value 5.5) • Alkalis/strong alkalis (neutral/enzymatic (max. permissible pH value 8.5, mandatory for devices made of aluminium or other alkali-sensitive materials) or alkaline cleaning agents (max. permissible pH value 11, mandatory for devices intended for use in prion-critical areas, e.g. in accordance with Appendix 7 of the KRINKO RKI BfArM recommendation for reprocessing)) • Organic solvents (e.g. alcohols, ethers, ketones, benzines) • Oxidising agents (e.g. hydrogen peroxides) • Halogens (chlorine, iodine, bromine) • Aromatic/halogenated hydrocarbons

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