


AUTOPSIESÄGE „ELECTRONIC POWER“  
AUTOPSY SAW „ELECTRONIC POWER“  
SIERRA DE AUTOPSIA „ELECTRONIC POWER“  
SCIE D'AUTOPSIE „ELECTRONIC POWER“



Grau dargestellte Artikelnummern bezeichnen „distributed by ELCON“ Produkte  
Article numbers shown in grey indicate „distributed by ELCON“ products  
Los numeros de articulos mostrados en gris indican productos „distribuidos por ELCON“  
Les numéros d'articles affichés en gris indiquent les produits „distribués par ELCON“

## AUTOPSY SAW, OSCILLATING



**550-004** Autopsy Saw, **standard** connector, 230 V / 50-60 Hz  
**550-006** Autopsy Saw, **standard** connector, 110 V / 50-60 Hz

complete with:

- 1 x sawblade Ø 50 mm
- 1 x sawblade Ø 65 mm
- 1 x sawblade Ø 65 mm
- 1 x pair of wrenches 11 cm / 1 x T-handle hexagon screw driver / 1 x grease / 1 x instruction manual

### TECHNICAL DATA

- 500 Watt power
- 3500 to 24000 oscillations, infinitely variable
- Fully automatic power-compensation
- Slim motor design with silent motor 52 – 77 dB(A)
- Latest electronic system with various protection and security functions
- Soft Start at switch-on to minimize accident risk and to save all mechanical and electrical components
- Special positioning of the motor
- Robust construction
- Safety ON / OFF switch
- Weight only 1,25 kg without cable
- Vibration- and maintenance-free
- 5 meter cable with EU-connection, 230 Volt  
3 meter cable with US-connection, 110 Volt
- Double protection isolation according to class II, medical use of type BF
- Radio and TV suppressed
- EN 60601-1, EN 60601-1-2, IEC 60601-1, IEC 60601-1-2 inspected
- 2 years warranty / 1 year guarantee on engine and transmission part (in case of proper use)



**550-000** Autopsy Saw, **triangle** connector, 230 V / 50-60 Hz  
**550-002** Autopsy Saw, **triangle** connector, 110 V / 50-60 Hz

complete with:

- 1 x shaft to fit in triangle connector (550-025)
- 1 x sawblade Ø 50 mm
- 1 x sawblade Ø 65 mm
- 1 x sawblade Ø 65 mm
- 1 x pair of wrenches 14 cm / 1 x T-handle hexagon screw driver / 1 x grease / 1 x instruction manual

### TECHNICAL DATA

- 500 Watt power
- 3500 to 24000 oscillations, infinitely variable
- Fully automatic power-compensation
- Slim motor design with silent motor 52 – 77 dB(A)
- Latest electronic system with various protection and security functions
- Soft Start at switch-on to minimize accident risk and to save all mechanical and electrical components
- Special positioning of the motor
- Robust construction
- Safety ON / OFF switch
- Weight only 1,25 kg without cable
- Vibration- and maintenance-free
- 5 meter cable with EU-connection, 230 Volt  
3 meter cable with US-connection, 110 Volt
- Double protection isolation according to class II, medical use of type BF
- Radio and TV suppressed
- EN 60601-1, EN 60601-1-2, IEC 60601-1, IEC 60601-1-2 inspected
- 2 years warranty / 1 year guarantee on engine and transmission part (in case of proper use)


**SAW BLADES**, to work with standard saw 550-004 + 006 and triangle connector saw 550-000 + 002 (additional shaft needed)



04-MST-00010    Ø 45 mm  
04-MST-00004    Ø 50 mm  
04-MST-00007    Ø 65 mm



04-MST-00014    Ø 65 mm




04-MST-00016    Ø 65 mm



04-MST-00018    8 mm



04-MST-00020    20 mm



**Eccentric disc**, to limit cutting depth  
04-MST-00022

**SHAFT**, needed to work with triangle connector saw 550-000 and 550-002



550-025    Shaft, 39 mm



550-026    Shaft, 58 mm



550-027    Shaft, double

**SAW BLADES WITH SHAFT**, pre-assembled for 550-000 and 550-002



**Saw Blade**, incl. shaft 550-025

550-030    Ø 50 mm  
550-031    Ø 65 mm



**Saw Blade**, incl. shaft 550-025

550-032    Ø 65 mm




**Saw Blade**, incl. shaft 550-025

550-033    Ø 65 mm



**Saw Blade**, incl. shaft 550-025

550-034    8 mm  
550-035    20 mm



**Saw Blades** 400-462-1 / 550-040  
incl. shaft 550-027

550-036    Ø 45 mm



**Eccentric disc**, to limit cutting depth  
**Saw Blade**, incl. shaft 550-026

550-037    Ø 65 mm

# **ELCON**

*MEDICAL INSTRUMENTS GmbH*

D - 78532 Tuttlingen  
Tel: +49 7461 9281-0  
info@elcon-medical.de

Dr. Karl-Storz-Str. 26  
Fax: +49 7461 9281-29  
www.elcon-medical.com