


3
4
3



ELCON
MEDICAL INSTRUMENTS GmbH

CLR RETRAKTOR SYSTEM
CLR RETRACTOR SYSTEM
CLR RÉTRACTEUR SYSTEM
CLR RETRACTOR SYSTEM

- Lumbar retractor system
- Side load for optimal exposure
- Blades out of titanium, color-coded, radiolucent
- Large variety of valves and hooks
- Made in Germany



Die aktuellen Informationen bezüglich der Risikoklasse der ELCON Produkte finden Sie auf unserer Homepage unter „Produkt finden“

You can find the latest information on risk classes of ELCON products on our homepage under “Find product.”

Las informaciones actualizadas sobre la clase de riesgo de los productos ELCON se encuentran en nuestra página web bajo „Buscar producto“

Vous trouverez les informations actuelles concernant la classe de risque des produit ELCON sur notre page d'accueil, sous la rubrique „Trouver un produit“

RETRACTOR / HANDLES



CLR Retraktor Side-load CLR retractor side-load	item no.
	343-3151 *

CCR/CLR Handgriff Side-load mit schwarzem Radelgriff CCR/CLR handle for blades with black PPSU-handle	item no.
	343-3120 *



CCR/CLR Handgriff Side-load, kanuliert mit Luer-Lock, mit schwarzem Radelgriff CCR/CLR handle for blades, cannulated with Luer-Lock, with black PPSU-handle	item no.
	343-3120LL *



CCR/ CLR Tiefenmeßstab, natur CCR/CLR depth gauge instrument, natural	item no.
	343-3125 *



BLADES

Alle Valven aus Titanium, Farbe: blau, röntgendurchlässig, mit Side-load-Mechanismus

All blades out of Titanium, color: blue, radiolucent, with side-load mechanism



CLR-Valve, side load, 20mm breit , blau CLR blade, side load, width: 20mm, blue		
width	depth	item no.
20	35	343-3160-35T *
20	40	343-3160-40T *
20	45	343-3160-45T *
20	55	343-3160-55T *
20	65	343-3160-65T *



CLR-Valve, side load, 40mm breit , blau CLR blade, side load, width: 40mm, blue		
width	depth	item no.
40	35	343-3161-35T *
40	40	343-3161-40T *
40	45	343-3161-45T *
40	55	343-3161-55T *
40	65	343-3161-65T *
40	75	343-3161-75T *



CLR-Valve, side load, 60mm breit , blau CLR blade, side load, width: 60mm, blue		
width	depth	item no.
60	35	343-3162-35T *
60	40	343-3162-40T *
60	45	343-3162-45T *
60	55	343-3162-55T *
60	65	343-3162-65T *
60	75	343-3162-75T *
60	85	343-3162-85T *




Hemi-laminectomy Haken, blau Hemi-laminectomy hook, blue	
	item no.
klein - small	343-3163-01T *
groß - large	343-3163-02T *

BASKETS

CLR-Siebkorb aus Edelstahl mit Deckel CLR-Tray, stainless steal with lid	item no.
Maße / dimensions: 485 mm x 253 mm x 70 mm	561-0526-70*



CLR Set Side-load bestehend aus Produkten mit (*) gekennzeichnet Valven paarweise CLR Set side load consisting of products marked (*) Blades, 2 pcs each	item no.
	343-3150-01



ELCON

MEDICAL INSTRUMENTS GmbH

D - 78532 Tuttlingen
Tel: +49 7461 9281-0
info@elcon-medical.de

Dr. Karl-Storz-Str. 26
Fax: +49 7461 9281-29
www.elcon-medical.com