


3  
4  
3



**ELCON**  
MEDICAL INSTRUMENTS GmbH

CCR RETRAKTOR SYSTEM  
CCR RETRACTOR SYSTEM  
CCR RÉTRACTEUR SYSTEM  
CCR RETRACTOR SYSTEM

- Cervical retraction system
- Titanium color-coded, radiolucent
- Side load for optimal exposure
- Large variety of blades
- Made in Germany



Die aktuellen Informationen bezüglich der Risikoklasse der ELCON Produkte finden Sie auf unserer Homepage unter „Produkt finden“  
You can find the latest information on risk classes of ELCON products on our homepage under “Find product.”

Las informaciones actualizadas sobre la clase de riesgo de los productos ELCON se encuentran en nuestra página web bajo „Buscar producto“

Vous trouverez les informations actuelles concernant la classe de risque des produit ELCON sur notre page d'accueil, sous la rubrique „Trouver un produit“


**RETRACTORS**



CCR Basisretractor Side-Load mit Doppelgelenk, Endgelenk aus Titanium CCR retractor Side-Load, double hinge, distal hinge out of titanium	item no.
	<b>343-3101 *</b>



CCR Gegensperrer Side-load, Endgelenk aus Titanium CCR counter retractor Side-load, distal hinge out of titanium	item no.
	<b>343-3102 *</b>



CCR Basis-Gegensperrer Side Load, mit Doppelgelenk, Endgelenk aus Titanium CCR basic counter retractor, double hinge, distal hinge out of titanium	item no.
	<b>343-3103 *</b>

### HANDLES

CCR/CLR Handgriff Side-load mit schwarzem Radelgriff CCR/CLR handle for blades with black PPSU-handle	item no.
	<b>343-3120 *</b>



CCR/CLR Handgriff Side-load, kanuliert mit Luer-Lock, mit schwarzem Radelgriff CCR/CLR handle for blades, cannulated with Luer-Lock, with black PPSU-handle	item no.
	<b>343-3120LL *</b>



CCR/ CLR Tiefenmeßstab, natur CCR/CLR depth gauge instrument, natural	item no.
	<b>343-3125 *</b>



**BLADES**

Alle Valven aus Titanium, röntgendurchlässig mit Side-load-Mechanismus  
all blades out of Titanium, radiolucent, with side-load-mechanism



CCR-Valve 4-zahnig grün CCR-blade green with 4 teeth		
width	depth	item no.
19	25	<b>343-3110-25T</b>
19	30	<b>343-3110-30T</b>
19	35	<b>343-3110-35T</b>
19	40	<b>343-3110-40T</b>
19	45	<b>343-3110-45T</b>
19	50	<b>343-3110-50T</b>
19	55	<b>343-3110-55T</b>
19	60	<b>343-3110-60T</b>
19	65	<b>343-3110-65T</b>
19	70	<b>343-3110-70T</b>
19	75	<b>343-3110-75T</b>



CCR-Valve, stumpf, longitudinal CCR-blade, blunt, longitudinal		
width	depth	item no.
19	25	<b>343-3111-25T</b>
19	30	<b>343-3111-30T</b>
19	35	<b>343-3111-35T</b>
19	40	<b>343-3111-40T</b>
19	45	<b>343-3111-45T</b>
19	50	<b>343-3111-50T</b>
19	55	<b>343-3111-55T</b>
19	60	<b>343-3111-60T</b>
19	65	<b>343-3111-65T</b>
19	70	<b>343-3111-70T</b>
19	75	<b>343-3111-75T</b>



CCR-Valve, stumpf, longitudinal CCR-blade, blunt, longitudinal		
width	depth	item no.
24	25	<b>343-3115-25T</b>
24	30	<b>343-3115-30T *</b>
24	35	<b>343-3115-35T *</b>
24	40	<b>343-3115-40T *</b>
24	45	<b>343-3115-45T *</b>
24	50	<b>343-3115-50T *</b>
24	55	<b>343-3115-55T *</b>
24	60	<b>343-3115-60T *</b>
24	65	<b>343-3115-65T *</b>
24	70	<b>343-3115-70T *</b>
24	75	<b>343-3115-75T</b>




CCR Valve, 5-zahnig medial CCR Valve, 5 teeth, medial		
width	depth	item no.
24	25	<b>343-3116-25T</b>
24	30	<b>343-3116-30T *</b>
24	35	<b>343-3116-35T *</b>
24	40	<b>343-3116-40T *</b>
24	45	<b>343-3116-45T *</b>
24	50	<b>343-3116-50T *</b>
24	55	<b>343-3116-55T *</b>
24	60	<b>343-3116-60T *</b>
24	65	<b>343-3116-65T *</b>
24	70	<b>343-3116-70T *</b>
24	75	<b>343-3116-75T</b>




CCR Valve, stumpf, longitudinal CCR Valve, blunt longitudinal		
width	depth	item no.
24	25	<b>343-3117-25T</b>
24	30	<b>343-3117-30T</b>
24	35	<b>343-3117-35T *</b>
24	40	<b>343-3117-40T *</b>
24	45	<b>343-3117-45T *</b>
24	50	<b>343-3117-50T *</b>
24	55	<b>343-3117-55T *</b>
24	60	<b>343-3117-60T *</b>
24	65	<b>343-3117-65T *</b>
24	70	<b>343-3117-70T *</b>
24	75	<b>343-3117-75T *</b>

### BASKETS

CCR-Siebkorb aus Edelstahl mit Deckel   CCR-Tray, stainless steel with lid	item no.
Maße / dimensions: 485 mm x 253 mm x 30 mm	<b>561-0525-30 *</b>



CCR Set bestehend aus Produkten welche mit „*“ gekennzeichnet (je 1 Stück) CCR-Set consisting of products marked with a „*“ (one piece each)	item no.
	<b>343-3100-01</b>



# **ELCON**

*MEDICAL INSTRUMENTS GmbH*

D - 78532 Tuttlingen  
Tel: +49 7461 9281-0  
info@elcon-medical.de

Dr. Karl-Storz-Str. 26  
Fax: +49 7461 9281-29  
www.elcon-medical.com