


3
4
3




ELCON
MEDICAL INSTRUMENTS GmbH

MCC RETRAKTOR
MCC RETRACTOR
MCC RÉTRACTEUR
MCC RETRACTOR

Micro Discectomy Retractor

- Optimal Exposure
- Retractors and blades TiAIN coated with black finish
- Made in Germany



Die aktuellen Informationen bezüglich der Risikoklasse der ELCON Produkte finden Sie auf unserer Homepage unter „Produkt finden“
You can find the latest information on risk classes of ELCON products on our homepage under “Find product.”

Las informaciones actualizadas sobre la clase de riesgo de los productos ELCON se encuentran en nuestra página web bajo „Buscar producto“

Vous trouverez les informations actuelles concernant la classe de risque des produit ELCON sur notre page d'accueil, sous la rubrique „Trouver un produit“

RETRAKTOREN, VALVEN: Titan beschichtet

RETRACTORS, BLADES: Titan coated



MCC Retraktor mit Gelenk, max. Spannweite 85 mm MCC Hinged Retractor, max. opening: 85 mm / 3.3"	item no.
	343-611 * / **



MCC Retraktor starr, max. Spannweite 85 mm MCC Rigid Retractor, max. opening: 85 mm / 3.3"	item no.
	343-612



Valve 2-zahnig, 17 mm breit Blade 2 teeth, 17 mm		
width	depth	item no.
17	30	343-617-30
17	40	343-617-40
17	50	343-617-50
17	60	343-617-60
17	70	343-617-70
17	80	343-617-80



Valve verjüngt, 6-zahnig, 27 mm breit Tapered blade, 6 teeth, 27 mm		
width	depth	item no.
27	40	343-618-40

VALVEN, HAKEN: Titan beschichtet

BLADES, HOOKS: Titan coated



Valve 4-zahnig, 20 mm breit Blade, 4 teeth, 20 mm		
width	depth	item no.
20	30	343-613-30 **
20	40	343-613-40 * / **
20	50	343-613-50 * / **
20	60	343-613-60 * / **
20	70	343-613-70 * / **
20	80	343-613-80 **



Valve 4-zahnig, 20 mm breit Blade, 4 teeth, 20 mm		
width	depth	item no.
27	30	343-614-30 **
27	40	343-614-40 * / **
27	50	343-614-50 * / **
27	60	343-614-60 * / **
27	70	343-614-70 * / **
27	80	343-614-80 **



Mehrzahn-Haken, 2 mm Multiple toothed hook , 2.0 mm		
Länge / length	Zähne / tooth	item no.
20	2	343-616-20 **
30	3	343-616-30 * / **
40	4	343-616-40 * / **
50	5	343-616-50 * / **
60	6	343-616-60 * / **
70	7	343-616-70 **



Haken, 1-zahnig, 1,8 mm Hook blade, 1 teeth, 1.8 mm		
Länge / length		item no.
20		343-615-20 **
30		343-615-30 * / **
40		343-615-40 * / **
50		343-615-50 * / **
60		343-615-60 * / **
70		343-615-70 **

SIEBKÖRBE / SETS

BASKETS / SETS


MCC Siebkorb klein mit Deckel MCC Tray small with lid	item no.
Maße / dimensions: 243 mm x 253 mm x 50 mm	561-0520-50 *

MCC Siebkorb groß mit Deckel MCC Tray large with lid	item no.
Maße / dimensions: 405 mm x 253 mm x 50 mm	561-0521-50 **

MCC Set klein bestehend aus den mit „ * “ gekennzeichneten Artikeln Valven paarweise, Haken einzeln Basic MCC Set consisting of the articles marked with „ * “ Blades 2 each, hook 1 each	item no.
	343-610-01



MCC Set groß bestehend aus den mit „ ** “ gekennzeichneten Artikeln Gelenkretraktor und Valven paarweise, Haken einzeln MCC Set large consisting of the articles marked with „ ** “ Hinged Retractor and Blades 2 each, hook 1 each	item no.
	343-610-02



ELCON

MEDICAL INSTRUMENTS GmbH

D - 78532 Tuttlingen
Tel: +49 7461 9281-0
info@elcon-medical.de

Dr. Karl-Storz-Str. 26
Fax: +49 7461 9281-29
www.elcon-medical.com